

Procainamide labelling as part of a sensitive workflow for the identification and quantitation of potential glycan drug targets in parasite cell membrane, insect vector saliva and host tissue models of infection by HILIC-HPLC-ESI-MS/MS

Katherine Wongtrakul-Kish, Radoslaw P. Kozak, Daniel I. R. Spencer
Ludger Ltd, Culham Science Centre, Oxfordshire, UK



Katherine.Wongtrakul-Kish@ludger.com

Introduction

Changes in glycosylation have been associated with many states of health and disease and have the potential to provide diagnostic information and act as drug targets. In parasitological studies however, investigations into glycosylation changes are hampered by the availability of small sample amounts. This problem is particularly challenging in analyses that focus on characterisation of plasma membrane glyco-conjugates. These glycans offer the most obvious choice of glycan-based drug targets due to their position at the interface between a cell and its environment. In such workflows however the enrichment of plasma membranes from an already small starting amount further decreases the amount of glycoprotein available.


Currently, the most widely used glycan fluorescence label in the analysis of mammalian glycosylation systems is 2-aminobenzamide (2-AB) which is compatible for use in HILIC-HPLC-FLR coupled with ESI-MS/MS. However 2-AB suffers from relatively poor MS signal intensity when analysing small sample amounts. For the analysis of novel or non-mammalian parasite glycans, MS signal intensity is an important factor in acquiring good quality, data-rich MS/MS fragmentation information. Here we evaluate the use of procainamide as an alternative to 2-AB in HILIC-HPLC-ESI-MS/MS glycomic analyses.

Aims:



- To determine a suitable glycomics workflow coupling HPLC with ESI-MS/MS to combat the problem of small sample amounts. Specifically, to compare and evaluate the use of two different labels (**2-AB and procainamide**) for released *N*-glycans from a standard glycoprotein (bovine fetuin) to ensure highly sensitive and reliable MS and FLR detection
- To investigate the role of *N*-glycosylation in various parasite, vector and host models of parasitological infection

Methods & Results


1. Comparison of 2-AB and procainamide labelled bovine fetuin *N*-glycans



2. Application of Procainamide labelling with HILIC-HPLC-FLR and ESI-MS/MS to the study of three model systems



3. Improved MS signal intensity provides high quality MS/MS fragmentation which contains diagnostic ions to aid in *N*-glycan structural characterisation



Conclusions

- Here we present a procainamide labelling workflow suitable for HPLC-MS analyses
- The workflow offers a sensitive and robust way of characterising and quantifying mammalian and non-mammalian glycans from small starting sample amounts.
- This workflow can be applied to a wide range of studies and biological samples.

Acknowledgements

This work was supported by funding from the GlycoPar EU FP7 Marie Curie Initial Training Network. We thank the following Network members and affiliates for the samples described here: Dr Sandra Trindade, Dr Joao Rodrigues, Dr Luisa Figueiredo, Dr Alvaro Acosta-Serrano, Karina Mondragon-Shem

References

- http://www.raywilsonbirdphotography.co.uk/Galleries/Invertebrates/vectors/Tsetse_Fly.html
- <http://www.taconic.com/mouse-model/black-6-b6ntac>
- <http://www.york.ac.uk/cii/media-library/>

